



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

Condition Report

Sale 1828 - 14/01/2019 10:00 - Live Webcast General Medical Equipment and Endoscopy Auction - PART 1

Sale Date 14-01-2019 10:00

Location: Malmesbury, Wiltshire, UK

Viewing: Friday 11th January 12PM - 5PM and morning of the sale from 08:30AM

Removal: By 4PM on Monday 21st January 2019 - Late collection charges will be incurred after this date

These reports are compiled by our Medical Engineering Team and are for guidance only - all lots are sold as seen, as per our Conditions of Sale. We recommend that wherever possible you inspect lots prior to purchase.

Lot 1 Olympus CF-P20S Fiber Sigmoidoscope | SN 2615735 | AN 358254

- Condition Report:
- **Insertion Tube** – Satisfactory condition
 - **Bending mechanism & Angulation Control** – Good movements and locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good
 - **Remote Control Switches** – N/A
 - **Distal End** – Good condition (distal cap missing)
 - **Image Quality** – Good image (with signs of a few broken fibres)
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Overall Condition** – Satisfactory overall physical condition with signs of wear on the handle
 - **Dimension:** 12 x 600 mm

Lot 2 Olympus GIF-2T240 Video Gastroscope | SN 2531261

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements and locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – Functional
 - **Overall Condition** – Good overall condition
 - **Dimension:** 12 x 1000 mm

Lot 3 Olympus CF-P20S Fiber Sigmoidoscope | SN 2536601



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

- Condition Report:
- **Insertion Tube** – Satisfactory condition
 - **Bending mechanism & Angulation Control** – Good movements and locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good
 - **Remote Control Switches** – N/A
 - **Distal End** – Good condition
 - **Image Quality** – Good image (with signs of a few broken fibres)
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Overall Condition** – Satisfactory overall condition with signs of wear on the handle
 - **Dimension:** 12 x 600 mm
-

Lot 4 Olympus CF-100S Video Sigmoidoscope | SN 2220762 - Item in Pentax Case

- Condition Report:
- **Insertion Tube** – Satisfactory condition
 - **Bending mechanism & Angulation Control** – Good movements (R/L lock is slightly loose)
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Not tested
 - **Remote Control Switches** – Not tested
 - **Distal End** – Satisfactory condition
 - **Image Quality** – Not tested
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Fail
 - **Overall Condition** – Satisfactory physical condition but leaking. Processor not available to check image quality. Auxiliary water connector is damaged.
 - **Dimension:** 13 x 600 mm
-

Lot 5 Fujinon FO-120S Choledocho Fibre Optic Scope | AN 038087 | SN HY124A002

- Condition Report:
- **Insertion Tube** – Bad kink on the edge of the control box
 - **Bending mechanism & Angulation Control** – Good movements
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good with 8 damaged fibres
 - **Remote Control Switches** – N/A
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good image but with kink on the insertion tube
 - **Dimension:** 6.2 x 400 mm
-

Lot 6 Olympus GIF-XQ260 Video Gastroscope | SN 2723491

- Condition Report:
- **Insertion Tube** – Good condition



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

- **Bending mechanism & Angulation Control** – Limited down angulation, locks ok
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition
 - **Dimension:** 10 x 1000 mm
-

Lot 7 Olympus CF-Q260 DL Video Colonoscope | SN 2710625

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition
 - **Dimension:** 12.5 x 1600 mm
-

Lot 8 Olympus CF-Q260 DL Video Colonoscope | SN 2710633

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition
 - **Dimension:** 12.5 x 1600 mm
-

Lot 9 Olympus CF-240DL Video Colonoscope | SN 2200060 (No Case)

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

- **Objective Lens & Focussing** – Good
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition
 - **Dimension:** 14.3 x 1600 mm
- No Carry case

Lot 10 Storz 13211 Fiberscope (Scope Type Unknown, No Case)

- Condition Report:
- **Insertion Tube** – Average condition
 - **Bending mechanism & Angulation Control** – Fair movements & locks
 - **Light Guide Cable & Light Transmission** – Fair
 - **Objective Lens & Focussing** – Poor optics
 - **Remote Control Switches** – N/A
 - **Distal End** – Fair condition
 - **Image Quality** – No image visible, poor
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – N/A
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Poor image and average physical condition
 - **Dimension:** 12 x 900 mm

Lot 11 Olympus GIF 2T100 Video Gastroscope | SN 2240253 (No Case)

- Condition Report:
- **Insertion Tube** – Fair condition
 - **Bending mechanism & Angulation Control** – limited down movements, locks ok
 - **Light Guide Cable & Light Transmission** – Fair
 - **Objective Lens & Focussing** – Not tested
 - **Remote Control Switches** – Not tested
 - **Distal End** – Good condition
 - **Image Quality** – Not tested
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Failed
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Leaking and in and average physical condition, image quality not tested
 - **Dimension:** 14 x 1000 mm
- No Carry case

Lot 12 Olympus CF 1T100L Video Colonoscope | SN 2552738 (No Case)

- Condition Report:
- **Insertion Tube** – De-laminated
 - **Bending mechanism & Angulation Control** – Fair movements & locks



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

- **Light Guide Cable & Light Transmission** – Fair
 - **Objective Lens & Focussing** – Not tested
 - **Remote Control Switches** – Not tested
 - **Distal End** – Fair condition
 - **Image Quality** – Not tested
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Failed
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Leaking and in and average physical condition, image quality not tested
 - **Dimension:** 14.6 x 1600 mm
- No Carry case

Lot 13 Olympus GIF 1T100 Video Gastroscope | SN 2100138 (No Case)

- Condition Report:
- **Insertion Tube** – Fair condition
 - **Bending mechanism & Angulation Control** – limited down movements, locks ok
 - **Light Guide Cable & Light Transmission** – Fair
 - **Objective Lens & Focussing** – Not tested
 - **Remote Control Switches** – Not tested
 - **Distal End** – Good condition
 - **Image Quality** – Not tested
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Average physical condition, image quality not tested
 - **Dimension:** 12.4 x 1000 mm
- No Carry case

Lot 14 Olympus CF-240DL Video Colonoscope | SN 2200059 (No Case)

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Failed
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Leaking but good physical condition and good image quality
 - **Dimension:** 12.5 x 1600 mm
- No Carry case

Lot 15 Olympus CF-1T20L Fiber Colonoscope | SN 2264153 (No Case)



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

- Condition Report:
- **Insertion Tube** – De-laminated
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Fair optics with 4 damaged fibres
 - **Remote Control Switches** – N/A
 - **Distal End** – Fair condition
 - **Image Quality** – Fair image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Fair overall condition
 - **Dimension:** 14.4 x 1600 mm
- No Carry case

Lot 16 Olympus Endoeve 30□ Video Scope

Condition Report: The instrument was tested on a compatible stack system, good light transmission and very good image quality displayed. Lot is in good physical condition. Tray WA59380A included.

Lot 17 DP Medical TFG Nasendoscope and Leakage Tester

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements
 - **Light Guide Cable & Light Transmission** – Good, fibre optic guide cable not included
 - **Objective Lens & Focussing** – Good opticsd
 - **Remote Control Switches** – N/A
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – N/A
 - **Leak Test** – Failed
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good image but leaking
 - **Dimension:** 4 x 300 mm

Lot 18 Machida LY-C30 Rigid Laryngoscope and Light Lead | SN 61028

Condition Report: The fibre optic guide cable is in good condition and the telescope's optics is clear with good magnification.

Lot 19 Machida LY-C30 Rigid Laryngoscope and Light Lead | SN 60533

Condition Report: The fibre optic guide cable is in good condition and the telescope's optics is clear with good magnification.

Lot 20 Machida LY-CS30 Rigid Laryngoscope and Light Lead | SN 60926



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

Condition Report: The fibre optic guide cable is in good condition and the telescope's optics is clear with good magnification.

Lot 21 Olympus CYF-240 Video Cystoscope | SN 2900616 (No Case)

Condition Report:

- **Insertion Tube** – Good condition
- **Bending mechanism & Angulation Control** – Good movements & lock
- **Light Guide Cable & Light Transmission** – Good
- **Objective Lens & Focussing** – Good
- **Remote Control Switches** – Fully functional
- **Distal End** – Good condition
- **Image Quality** – Good image
- **Instrument Channel (s)** – Clear
- **Leak Test** – Pass
- **Forceps Elevation Mechanism** – N/A
- **Overall Condition** – Good overall condition
- **Dimension:** 5.6 x 350 mm

No Carry case

Lot 22 Olympus CYF-240 Video Cystoscope | SN2800610 (No Case)

Condition Report:

- **Insertion Tube** – Good condition
- **Bending mechanism & Angulation Control** – Good movements & lock
- **Light Guide Cable & Light Transmission** – Good
- **Objective Lens & Focussing** – Good
- **Remote Control Switches** – Fully functional
- **Distal End** – Good condition
- **Image Quality** – Good image
- **Instrument Channel (s)** – Clear
- **Leak Test** – Pass
- **Forceps Elevation Mechanism** – N/A
- **Overall Condition** – Good overall condition
- **Dimension:** 5.6 x 350 mm

No Carry case

Lot 23 Olympus CYF-240 Video Cystoscope | SN 2800614 (No Case)

Condition Report:

- **Insertion Tube** – Good condition
- **Bending mechanism & Angulation Control** – Good movements & lock
- **Light Guide Cable & Light Transmission** – Good
- **Objective Lens & Focussing** – Good
- **Remote Control Switches** – Fully functional
- **Distal End** – Good condition
- **Image Quality** – Good image
- **Instrument Channel (s)** – Clear
- **Leak Test** – Pass



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

- **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition
 - **Dimension:** 5.6 x 350 mm
- No Carry case

Lot 24 Olympus CYF-240 Video Cystoscope | SN 2900696 (No Case)

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & lock
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition
 - **Dimension:** 5.6 x 350 mm
- No Carry case

Lot 25 Olympus ENF XP Fibre Rhinolaryngoscope | SN 1100665 (No Case)

- Condition Report:
- **Insertion Tube** – Crush on the bending section
 - **Bending mechanism & Angulation Control** – Good movements
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Poor optics with several damaged fibres
 - **Remote Control Switches** – N/A
 - **Distal End** – Fair condition
 - **Image Quality** – Poor image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Failed
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Poor image with damaged fibres, leaking and crush on the bending section
 - **Dimension:** 2 x 300 mm
- No Carry case

Lot 26 Olympus CF-Q260DL Video Colonoscope | SN 2411932

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good image



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

- **Instrument Channel (s)** – Clear
- **Leak Test** – Pass
- **Forceps Elevation Mechanism** – N/A
- **Overall Condition** – Good overall condition
- **Dimension:** 12.5 x 1600 mm

Lot 27 Olympus CF-240DL Video Colonoscope | SN 2200068

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Failed
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good image quality but leaking
 - **Dimension:** 12.5 x 1600 mm

Lot 28 Olympus TJF-200 Video Duodenoscope | SN 2400386

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Fair optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Fair image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – Functional
 - **Overall Condition** – Good overall condition
 - **Dimension:** 12 x 1200 mm

Lot 29 Olympus CF-Q260DL Video Colonoscope | SN 2211541

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

- **Forceps Elevation Mechanism** – N/A
- **Overall Condition** – Good overall condition
- **Dimension:** 14 x 1600 mm

Lot 30 Olympus CF-Q260DL Video Colonoscope | SN 2610475

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition
 - **Dimension:** 14 x 1600 mm

Lot 31 Olympus BF-UC 260FW Video Bronchoscope | SN 1010208

- Condition Report:
- **Insertion Tube** – Fair condition
 - **Bending mechanism & Angulation Control** – Good movements
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Fair optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Average image quality
 - **Instrument Channel (s)** – N/A
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Reasonable image quality. Ultrasound function was not tested
 - **Dimension:** 6.4 x 600 mm

Lot 32 Olympus CF-H260AZL Video Colonoscope | SN 2800184

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

- **Dimension:** 14.2 x 1600 mm
-

Lot 33 Olympus CF-H260AZL Video Colonoscope | SN 2700142

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Very good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Failed
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Very good image quality but with small leak
 - **Dimension:** 14 x 1600 mm
-

Lot 34 Olympus TJF-240 Video Duodenoscope | SN 2030298

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Fair optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition, end cap not fitted
 - **Image Quality** – Fair image quality
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Failed
 - **Forceps Elevation Mechanism** – Functional
 - **Overall Condition** – Fair image quality but leaking
 - **Dimension:** 12 x 1200 mm
-

Lot 35 Olympus GIF-1T240 Video Gastroscope | SN 2030051

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Limited down movement, locks ok
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition
 - **Dimension:** 12.2 x 1000 mm
-



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

Lot 36 Olympus BF-240 Video Bronchoscope | SN 1040453

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition
 - **Dimension:** 6.2 x 500 mm
-

Lot 37 Olympus GIF-Q240Z Video Gastroscope | SN 2000095

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Fair movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Average optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Image is a bit blurry
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Failed
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Blurry image and leaking
 - **Dimension:** 11.5 x 1000 mm
-

Lot 38 Olympus CF-Q260DL Video Colonoscope | SN 2610470

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition with scratches on control handle
 - **Dimension:** 14 x 1600 mm
-

Lot 39 Olympus GIF-H260 Video Gastroscope | SN 2948580



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Very good Image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition with scratches on control handle
 - **Dimension:** 11 x 1000 mm
-

Lot 40 Olympus GIF-2T240 Video Gastroscope | SN 2931028

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Fair movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good Image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition with scratches on control handle
 - **Dimension:** 13 x 1000 mm
-

Lot 41 Olympus GIF-H260Z Video Gastroscope | SN 2810814

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Fair movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Very good Image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition with scratches on control handle
 - **Dimension:** 12.2 x 1000 mm
-

Lot 42 Olympus GIF-IT240 Video Gastroscope | SN 2030055

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Fair movements & locks



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

- **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good Image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition with scratches on control handle
 - **Dimension:** 12.2 x 1000 mm
-

Lot 43 Olympus GIF-IT240 Video Gastroscope | SN 2030057

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Fair movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good Image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition with scratches on control handle
 - **Dimension:** 12.2 x 1000 mm
-

Lot 44 Olympus GIF-H260Z Video Gastroscope | SN 2911414

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Fair movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good Image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition
 - **Dimension:** 12.2 x 1000 mm
-

Lot 45 Olympus CF-H260AZL Video Colonoscope | SN 2900230

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

- **Remote Control Switches** – Fully functional
- **Distal End** – Good condition
- **Image Quality** – Good image
- **Instrument Channel (s)** – Clear
- **Leak Test** – Pass
- **Forceps Elevation Mechanism** – N/A
- **Overall Condition** – Good overall condition with minor scratches on control handle
- **Dimension:** 13.8 x 1600 mm

Lot 46 Olympus GIF-2TQ260M Video Gastroscope | SN 2810486

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good Image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition
 - **Dimension:** 13 x 1000 mm

Lot 47 Olympus TJF-240 Video Duodenoscope | SN 2431011

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Fair movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Fair condition, end cap not fitted
 - **Image Quality** – Good Image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Failed
 - **Forceps Elevation Mechanism** – Functional
 - **Overall Condition** – Good image but leaking
 - **Dimension:** 12.6 x 1200 mm

Lot 48 Olympus CF-240DL Video Colonoscope | SN 2200067

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

- **Image Quality** – Good image
- **Instrument Channel (s)** – Clear
- **Leak Test** – Pass
- **Forceps Elevation Mechanism** – N/A
- **Overall Condition** – Good overall condition with scratches on control handle
- **Dimension:** 14.2 x 1600 mm

Lot 49 Olympus CF-Q240DL Video Colonoscope | SN 2411931

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition with scratches on control handle
 - **Dimension:** 14 x 1600 mm

Lot 50 Olympus GIF-H260 Video Gastroscope | SN 2344946

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Fair movements & locks
 - **Light Guide Cable & Light Transmission** – Good
 - **Objective Lens & Focussing** – Fair optics
 - **Remote Control Switches** – Fully functional
 - **Distal End** – Good condition
 - **Image Quality** – Average Image quality, white balance failed
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Poor image, unable to adjust white balance
 - **Dimension:** 11 x 1000 mm

Lot 51 Olympus LF-GP Tracheal Intubation Fiberscope | SN W800626

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements
 - **Light Guide Cable & Light Transmission** – Good but fibre optic guide cable is not included
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – N/A
 - **Distal End** – Good condition
 - **Image Quality** – Good image



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

- **Instrument Channel (s)** – Clear
- **Leak Test** – Pass
- **Forceps Elevation Mechanism** – N/A
- **Overall Condition** – Good overall condition
- **Dimension:** 4.5 x 550 mm

Lot 52 Olympus LF-TP Tracheal Intubation Fiberscope | SN W500058

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements
 - **Light Guide Cable & Light Transmission** – Good but fibre optic guide cable is not included
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – N/A
 - **Distal End** – Good condition
 - **Image Quality** – Fair image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Pass
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good overall condition
 - **Dimension:** 4.5 x 550 mm

Lot 53 Olympus LF-GP Tracheal Intubation Fiberscope | SN 1414401

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements
 - **Light Guide Cable & Light Transmission** – Good but fibre optic guide cable is not included
 - **Objective Lens & Focussing** – Fair optics with 14 damaged fibres
 - **Remote Control Switches** – N/A
 - **Distal End** – Good condition
 - **Image Quality** – Average image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Failed
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Fair image but leaking
 - **Dimension:** 4.5 x 550 mm

Lot 54 Olympus LF-GP Tracheal Intubation Fiberscope | SN W700521

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements
 - **Light Guide Cable & Light Transmission** – Good but fibre optic guide cable is not included



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

- **Objective Lens & Focussing** – Good optics
- **Remote Control Switches** – N/A
- **Distal End** – Good condition
- **Image Quality** – Good image
- **Instrument Channel (s)** – Clear
- **Leak Test** – Pass
- **Forceps Elevation Mechanism** – N/A
- **Overall Condition** – Good overall condition
- **Dimension:** 4.5 x 550 mm

Lot 55 Olympus LF-GP Tracheal Intubation Fiberscope | 1414171

Serial No: 1414171

- Condition Report:
- **Insertion Tube** – Good condition
 - **Bending mechanism & Angulation Control** – Good movements
 - **Light Guide Cable & Light Transmission** – Good but fibre optic guide cable is not included
 - **Objective Lens & Focussing** – Good optics
 - **Remote Control Switches** – N/A
 - **Distal End** – Good condition
 - **Image Quality** – Good image
 - **Instrument Channel (s)** – Clear
 - **Leak Test** – Failed
 - **Forceps Elevation Mechanism** – N/A
 - **Overall Condition** – Good image quality but leaking
 - **Dimension:** 4.5 x 550 mm

Lot 73 1 x GE E-CAiov Gas Module and 1 x GE E-PRESTN Module with Quantity of Leads

Serial No: 6424019 (E-CAiOV) & 6473009 (E-PRESTN)

Manufactured: 2008

Accessories: 5-wire ECG lead x1
Masino SpO2 finger clip sensor
NIBP hose & cuff

Condition Report: The multiparameter module E-PRESTN was tested on a compatible S5 monitor and patient simulator. Good graphic and parameter values were displayed as per simulator set values. The Gas module failed to work when plugged in. Initial gas calibration start follow by "Sensor INOP" is displayed. Both modules are in fair physical condition. The SpO2 sensor included is not working.

Lot 74 1 x GE E-CAiov Gas Module and 1 x GE E-PRESTN Module with Quantity of Leads

Serial No: 6466927 (E-CAiOV) & 6472997 (E-PRESTN)

Manufactured: 2008

Accessories: 5-wire ECG lead x1
Masino SpO2 finger clip sensor



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

NIBP hose & cuff

Condition Report: The multiparameter module E-PRESTN was tested on a compatible S5 monitor and patient simulator. Good graphic and parameter values were displayed as per simulator set values. The Gas module was tested with calibration gas and all the readings were slightly below the expected values (need cal). Both modules are in fair physical condition. The SpO2 sensor included is not working.

Lot 186 Pentax EPK-1000 Processor | DOM 2010 | SN EDO14318

Manufactured: 2010

Condition Report: The unit powers up and the front panel controls and indicators are functional. The xenon lamp and brightness control are in good working order and the air pump is also working. The unit has not been fully tested with a video endoscope but the video output was present when connected to a monitor through an s-video cable. Overall the unit is in reasonable physical condition with minor signs of wear.

Lot 189 Valleylab Force E2 Diathermy | SN F1H5373B (No Foot Pedals)

Serial No: F1H5373B

Manufactured: 06 2001

Condition Report: Unit power up with the alphanumerical displays and function keys working fine. The display power up test completed. Patient return plate monitoring circuit is functional. The system was tested with a ESU analyser in both bipolar and monopolar functions and different modes. All the RF outputs were within the set limits. Foot switches are not included. The unit is in fair physical condition with some minor marks and scratches.

Lot 190 Coviden Valley Lab Force FXC Diathermy (No Foot Pedals) | DOM 07/2005 | SN

Serial No: F5G40135A

Manufactured: 07 2005

Condition Report: Unit power up with the alphanumerical displays and function keys working fine. The display power up test completed. Patient return plate monitoring circuit is functional. The system was tested with a ESU analyser in both bipolar and monopolar functions and different modes. All the RF outputs were within the set limits. Foot switches are not included. The unit is in good physical condition.

Lot 212 Fuji Fujinon 4400 System Processor | DOM 2007 | SN 7V396A156, 4400 System

Condition Report: The processor and light source both power up and the front panel controls and indicators are functional. The xenon lamp and brightness control are in good working order and the air pump is also working. Light source interface cable is included. The system has not been fully tested with a video endoscope but the video output was present when connected to a monitor through an s-video cable. Overall the units are in good physical condition with minor scratch marks on the light source.



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

Lot 213 Stryker 1188 High Def Camera Unit | DOM 03/2008 | SN 08C018384 and X3000

Condition Report: The processor and light source both power up and the front panel controls and display screens are functional. The xenon lamp and brightness control are in good working order. The system was tested with a compatible camera head and clear images were produced on the test monitor. Overall the units are in reasonable physical condition with minor signs of wear.

Lot 214 Olympus CV-260 Camera Processor | SN 7525133 and CLV-260 Light Source | SN

Condition Report: The processor and light source both power up and the front panel controls and indicators are functional. The xenon lamp and brightness control are in good working order and the air pump is also working. Light source interface cable is not included. The system has not been fully tested. Overall the units are in satisfactory physical condition with minor signs of wear. The front covers of the light source and processor are slightly discoloured.

Lot 218 Olympus Evis Lucera CV-260 Processor with MAJ-1154 Pigtail Video Connector | SN

Condition Report: The light source powers up and the front panel controls and indicators are functional. The xenon lamp and brightness control are in good working order and the air pump is also working. The processor powers up and the front panel indicators illuminate briefly but the system does not respond to any of the controls. The pigtail is in good working order (tested with a compatible endoscopy stack system and video scope). Light source interface cable is not included. Overall the units are in reasonable physical condition with minor signs of wear. The front cover of the processor is slightly discoloured. The pigtail has scratch marks and a small chip on one end of the connector.

Lot 264 Stryker 1088HD Camera Head | SN 06B005814

Condition Report: The camera head is in good working order, tested with a compatible control unit and the system produced clear images. Focus control and white balance setting are also functional. The image showed slight interference when the camera head cable was shaken. Otherwise the camera head is in reasonable physical condition.

Lot 265 Stryker 1088HD Camera Head | SN 06B003794

Condition Report: The camera head is in good working order, tested with a compatible control unit and the system produced clear images. Focus control and white balance setting are also functional. Overall the camera head is in reasonable physical condition.



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

Lot 266 Stryker 1188HD Camera | SN 07A005754

Condition Report: The camera head is in good working order, tested with a compatible control unit and the system produced clear images. Focus control and white balance setting are also functional. The image showed slight interference when the camera head cable was shaken. Otherwise the camera head is in reasonable physical condition.

Lot 267 Stryker 1188HD Camera | SN 10G014754

Condition Report: The camera head is in good working order, tested with a compatible control unit and the system produced clear images. Focus control and white balance setting are also functional but the buttons are missing from the camera head. The image showed slight interference when the camera head cable was shaken. Otherwise the camera head is in reasonable physical condition.

Lot 375 Storz Stack Trolley with SCB Electronic Endoflator | SN GC15071-B | Xenon Nova

Condition Report: All units within the stack system power up accordingly. The light source is fully functional with xenon lamp and brightness control in good working order. The Image 1 hub processor appears to be functional but has not been fully tested (no camera head available). The colour bar video output is present on the LCD monitor when the system is on. The insufflator, printer and DVD recorder all appear to be working but have not been fully tested. Overall the units are in reasonable physical condition. The stack trolley does not move smoothly (one of the wheels is damaged).

Lot 376 Olympus Stack Trolley with OEV261H Monitor | SN 7257783 | DOM 11/20121 | CV

Condition Report: The processor and light source both power up and the front panel controls and indicators are functional. The xenon lamp and brightness control are in good working order and the air pump is also working. The LCD monitor also powers up and appears to be working. The system has not been fully tested (some of the interconnecting cables including the light source interface cable are missing). Overall the units are in reasonable physical condition with minor signs of wear.

Lot 377 Stryker Stack Trolley with 2 x SV-2 High Def Monitors, X7000 Light source | SN

Condition Report: The light source is fully functional with xenon lamp and brightness control in good working order. The 1088i HD processor and camera head are also working and the



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

system produces clear images. White balance setting and focus control of the camera head are also functional. Both LCD monitors are working. The SDC HD image management unit powers up but displays diskboot failure (HDD appears to be missing). The printer has not been tested. Overall the units are in satisfactory physical condition. One of the monitors has scratch marks on the LCD protector.

Lot 378 Storz Urban M1057 Operating Microscope with Eyepieces and Lens | SN 0250

Condition Report: The operating microscope powers up, lamp and brightness control are in good working order. The optical system is clear, the binocular assembly is intact with eyepieces and objective lens in place and the system produces clear images. The focus and magnification controls are functional. The column height adjustment and lock are also working. Overall the unit is in reasonable physical condition with minor scratch marks on the control panel. The base locks are not holding firmly.
Objective lens = 250mm. Eyepieces = 10x.

Lot 379 Wild Heerbrugg M650 Operating Microscope with Lens and Eyepieces

Condition Report: The operating microscope powers up but the lamp is not working (needs replacing). The optical system appear to be clear. The binocular assembly is intact with eyepieces and objective lens in place. The focus and magnification controls are functional. Overall the unit is in reasonable physical condition with scratch marks on the paintwork.
Objective lens = 250mm. Eyepieces = 10x/21.

Lot 380 EUK Stack Trolley with Fujinon System 4400 Light Source | DOM 2008 | SN

Condition Report: The processor and light source both power up and the front panel controls and indicators are functional. The xenon lamp and brightness control are in good working order and the air pump is also working. The processor also appears to be functional but the system has not been fully tested (no endoscope available). The LCD monitor is working and the video output from the processor is present on the LCD monitor when the system is on. Overall the units are in satisfactory physical condition. The outer casing of the light source and processor have signs of corrosion and minor discolouration to the front covers. The rear panel of the stack trolley is not secure.

Lot 382 Datascope System 98 Balloon Pump with Datascope IABP Doppler

Operating Hours: 5158

Condition Report: The unit powers up and completes the initial system self test successfully. LCD display screen, control panel, pump and printer are in good working order. The



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

adult condensate removal / safety disk is in place (expiration date MArch 2013) but the Helium gas cylinder is not included. Doppler is in place but the battery appears to be depleted. No further checks carried out. Overall the unit is in reasonable physical condition.

Lot 384 Ethicon Mammotome Breast Biopsy System with Devicor Universal Targeting Set

Serial No: 2304195

Accessories: Devicor Universal Targeting Set - Long 145mm
Devicor MR Bladless Probe with Tube Set - Short 115mm
System Operating Manual CD-ROM
Mammotome Biopsy System Service Manual CD-ROM
Control Module Interlock

Condition Report: System powers up with the LCD working. System requires all devices to be connected. Unable to do any further tests.

Lot 426 Valley Lab Force FXC Diathermy (No Foot Pedals) | SN F3H28442A

Serial No: F3H28442A

Manufactured: 08 2003

Condition Report: Unit power up with the alphanumerical displays and function keys working fine. The display power up test completed. Patient return plate monitoring circuit is functional. The system was tested with a ESU analyser in both bipolar and monopolar functions and different modes. All the RF outputs were within the set limits. Foot switches are not included. The unit is in good physical condition.

Lot 427 Valley Lab Force FXC Diathermy (No Foot Pedals) | SN 2L24723A

Serial No: F2L24723A

Manufactured: 11 2002

Condition Report: Unit power up, but does not proceed through self-test. Bipolar display panel flashes and the unit emits an audible alarm. There are no footswitches included with the unit. The unit is in good overall physical condition.

Lot 435 Draeger Airshield Resuscitaire | DOM 2010 | DOM 2012 | SN YX06392 | SN

Accessories: Air & O2 hose

Condition Report: The unit powers up and the displays, function keys are working fine. The overhead heater and light are working fine when tested in the pre-warm and manual modes. The baby mode was not fully tested, skin temperature probe is not included. The motorised height adjustment is in good working order. The pneumatics system includes O2/air blender, AutoBreath, suction, Auxiliary flow & patient supply are in good working order. All flow meter and manometers are functional. The Cylinder pressure gauge has a small crack on the perspex but working fine. Mechanical cot tilt and mattress are in fair condition, with scratches on the handle. Overall the unit is in good physical condition.



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

Lot 438 Valley Lab Force FXC Diathermy (No Foot Pedals) | SN F4L35772A

Serial No: F4L35772A

Manufactured: 11 2004

Condition Report: Unit power up with the alphanumeric displays and function keys working fine. The display power up test completed. Patient return plate monitoring circuit is functional. The system was tested with a ESU analyser in both bipolar and monopolar functions and different modes. All the RF outputs were within the set limits. Foot switches are not included. The unit is in fair physical condition with some marks and slight damage to upwards key next to "Coag" display.

Lot 439 Valley Lab Force FXC Diathermy (No Foot Pedals) | SN F3L29939A

Serial No: F3L29939A

Manufactured: 11 2003

Condition Report: Unit power up with the alphanumeric displays and function keys working fine. The display power up test completed. Patient return plate monitoring circuit is functional. The system was tested with a ESU analyser in both bipolar and monopolar functions and different modes. All the RF outputs were within the set limits. Foot switches are not included. The unit is in fair physical condition with some marks and slight damage to upwards key next to "Coag" display.

Lot 481 Blease Frontline Genius Anaesthetic Machine with Blease 2200 MRI Ventilator

Serial No: 0840106

Associated Equipment: Blease 2200 MRI Ventilator, S/N: 0506-024

Accessories: Air, O2 & N2O hose

Condition Report: The anaesthetic system was tested with medical air and test lung. The flowmeters and pressure gauges are in good working order. The Blease 2200 ventilator is working fine with good ventilation cycles delivered. AlarmPAC alarm need battery replacement.. O2 monitor is fully functional. System is in good working order. Low O2 alarm is working. Lot is in good overall condition.

Lot 482 Draeger Evita 4 Ventilator Serial: ARRA-0135

Serial No: ARRA-0135

Accessories: Air & O2 hose

Software: Version 04.22

Condition Report: The Evita 4 power up and completed the initial power on tests. Touch screen, function keys and rotary selector knob working fine. The LCD screen is slightly darker at the bottom end. The function menus and ventilator modes are selectable. Exp. valve and flow sensor are not included but a set was used for testing the system. Good ventilation cycles were delivered with the Vt and frequency (BPM) corresponding to the measured minute volume (MV). Oxygen sensor was not working at the time of test, need cal. Vent mode available are: SIMV, BIPAP, ASB & APRV. Unit is in fair physical condition with minor scuffs. Trolley is not the original Draeger one. Backup internal battery is working.

Lot 515 Olympus OEV191H LCD Monitor, CV260SL Camera Unit, CLV-260SL Light



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

Condition Report: The processor and light source both power up and the front panel controls and indicators are functional. The xenon lamp and brightness control are in good working order and the air pump is also working. The LCD monitor, printer and DVD recorder also power up and appear to be working. The units are on a Fujinon stack trolley and the interconnecting cables are not included. The system has not been fully tested. Overall the units are in reasonable physical condition with minor signs of wear.

Lot 530 Olympus Stack Trolley with LCD Monitor, Evis Exera IU CV-180 Processor, CLV

Serial No: 7631981 (CV-180)

Associated Equipment: Sony Video Printer UP-21MD
Sony DVD Recorder DVO-100MD

Condition Report: Both CV-180 processor & CLV-180 light source power up with displays, indicators and key pad are working fine. The processor displayed good image quality when the camera head is connected. The light source's xenon lamp and brightness control are fully functional. LCD monitor and keyboard are working fine. The DVD recorder and printer are included but not fully tested. Lot is in good overall condition with minor scuffs and scratches consistent with use. Est. DOM: 2008.

Lot 531 Storz Stack System with Monitor, Aida Control | 20096020, Image 1 Hub |

Condition Report: The Image 1 Hub power up but was not fully tested, no camera head included. Electronic Endoflator power up with the display, flow 7 pressure led bars and key pad in good working order. System was not tested with CO2 cylinder. The xenon light source power up with lamp and brightness control fully functional. LCD with this lot. AIDA DVD unit power up but the CMOS battery needs replacing, does not retain the date settings.

Lot 532 Storz Stack System with Twin Monitors, SCB Image 1 Hub | 22200020, Aida DVD |

Condition Report: The Image 1 power up but was not fully tested, no camera head included. Electronic Endoflator power up with the display, flow 7 pressure led bars and key pad in good working order. System was not tested with CO2 cylinder. The xenon light source power up with lamp and brightness control fully functional. The main LCD monitor does not have mains adapter. There are a few scratches on the screen. AIDA DVD power up and touch screen is on but the DVD tray is stuck.

Lot 538 Olympus Stack Trolley with CV-260 Image Processor | SN 7647930, CLV-260 Light

Serial No: 7647930 (CV-260)

Manufactured: 2006

Associated Equipment: Sony Video Printer UP-55MD
Sony DVD Recorder DVO-1000MD
Olympus CLV-260 light source, S/N: 7617123

Accessories: Video connector MAJ1154
Inter link cable

Condition Report: All equipment powers up on both display trolley and instrument trolley. The video



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

processor indicator and function keys are working fine. System was tested with a compatible videoscope and good image was processed and displayed on the 19" LCD monitor. Light source's xenon lamp and brightness control are working fine. The video printer come on but it is out of ribbon. DVD recorder power but was not fully tested. System is in a fair physical condition.

Lot 565 Blease Genius Anaesthetic Machine with Datum Isoflurance Vapouriser and

Serial No: GENI-000003

Manufactured: 09 2008

Associated Equipment: Blease 2200 MRI ventilator, S/N: 20802-006

Accessories: Spacelab Rebreathing assembly
Blease Isoflurane vapouriser
AlarmPAC, model AP
Air, N2O & O2 hose
AGSS receiver

Condition Report: The anaesthetic system was tested with medical air and test lung. The flowmeters and pressure gauges are in good working order. The Blease 2200 ventilator settings and AlarmPAC alarm settings are working fine with good ventilation cycles delivered through the re-breathing assembly. O2 monitor is fully functional. System is in good working order. Low O2 alarm is working. Lot is in good overall condition.

Lot 566 Nikkiso DBB-05 Dialysis Machine | DOM 12/2009 | SN 69138N04 | AN 113602 and

Software: Version 2.0

Condition Report: Both units power up and completed the CPU & memory tests. Touch screen in good working order. The function menus are accessible but system was not fully tested. One of the units has a missing arterial blood pump rotor missing. Units are in fair physical condition. Operating hours are: 311188 & 210912.

Lot 567 Draeger Airshields Resuscitaire | DOM 2009 | SN TE05299

Manufactured: 2009

Accessories: Air & O2 hose

Condition Report: The unit powers up and the displays, function keys are working fine. The overhead heater and light are working fine when tested in the pre-warm and manual modes. The baby mode was not fully tested, skin temperature probe is not included. The motorised height adjustment is in good working order. The pneumatics system include O2/air blender, suction, Auxiliary flow & patient supply are in good working order. All flow meter and manometers are functional. The Cylinder pressure gauge has a small crack on the perspex but working fine. Mechanical cot tilt and mattress are in fair condition, with scratches on the handle. One of the perspex side cot (front) is missing, otherwise unit is in fair physical condition.

Lot 571 Blease Frontline Genius Anaesthetic Machine with Blease 2200 MRI Ventilator |

Accessories: Blease 2200 MRI Ventilator, S/N: 20009-002
AlarmPAC



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

Air, N2O & O2 hose
AGSS receiver

Condition Report: The anaesthetic system was tested with medical air and test lung. The flowmeters and pressure gauges are in good working order. The Blease 2200 ventilator settings and AlarmPAC alarm settings are working fine with good ventilation cycles delivered. O2 monitor is fully functional. System is in good working order. Low O2 alarm is working. Lot is in good overall condition.

Lot 573 Eschmann MR Hydraulic Operating Table | SN MR-1317 | AN 020940

Serial No: MR-1317

Condition Report: The mechanical adjustment for tilt, trendelenberg and reverse trendelenberg are working fine. Gas spring operated head rest and leg support are fully functional. Hydraulic lift is working fine. The table frame is in a fair condition with few scratches. Cushions are in fair condition. Casters and brake are working.

Lot 576 Datex Ohmeda S/5 Avance Anaesthetic Machine | DOM 2007 | SN ANBL01773 |

Accessories: Air, O2 & N2O hose

Condition Report: The S/5 Avance anaesthetic system power up and completed the electronic self-test. LCD, rotary selector knob and key pad are working. The system has an internal leak when air is connected and the bellow failed to pressurised. The work top cover is missing. The backup battery is not working.

Lot 578 Nikkiso DBB-05 Dialysis Machine | DOM 09/2005 | SN 65093-04 | AN 113342 and

Serial No: 65093-04 & 66158-02

Manufactured: 2005

Condition Report: Both units power up and completed the CPU & memory tests. Touch screen in good working order. The function menus are accessible but system was not fully tested. Units are in fair physical condition. Operating hours are: 336625 & 315730.

Lot 583 Draeger Airshield Resuscitaire | DOM 2009 | SN VE05539

Serial No: UE05539

Manufactured: 2009

Accessories: Air & O2 hose

Condition Report: The unit powers up and the displays, function keys are working fine. The overhead heater and light are working fine when tested in the pre-warm and manual modes. The baby mode was not fully tested, skin temperature probe is not included. The motorised height adjustment is in good working order. The pneumatics system include O2/air blender, AutoBreath, suction, Auxiliary flow & patient supply are in good working order. All flow meter and manometers are functional. The Cylinder pressure gauge has a small crack on the perspex but working fine. Mechanical cot tilt and mattress are in fair condition, with scratches on the handle. Overall the unit is in fair physical condition.

Lot 584 Nellcor Puritan Bennet Ventilator System | SN 3510080399

Serial No: 3510080399 (BDU)



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

Manufactured: 2008

Accessories: Air & O2 hose

Battery Pack S/N: 0816000407

Graphical User Interface Unit (GUI) S/N: 0817090461

Condition Report: The 840 ventilator system power up but the Graphic User Interface Unit does not come on. No further test done. Battery pack is working fine.

Lot 585 Nellcor Puritan Bennet Ventilator System | SN 3510041887

Serial No: 350041887 (BDU)

Manufactured: 2004

Operating Hours: 44375

Accessories: Air & O2 hose

Battery Pack S/N: 0416001923

Graphical User Interface Unit (GUI) S/N: 0417091917

Software: 4-070000-85-AN

Condition Report: The 840 ventilator system power up and completed the initial power on tests. GUI's LCD/touch screen key pad and rotary selector knob working fine. The function menus are accessible and the ventilator modes are selectable. Exp. filter is not included but one was used for testing the system. The BDU is working fine with good ventilation cycles were delivered. Tidal volume (VT) and frequency (BPM) corresponding to the measured minute volume (MV). Oxygen sensor is alarming but flow sensor is working. Vent mode available are: Spont, SIMV, Bilevel, A/C, PC, VC & PS. Software option installed: BiLevel, TC, VV+, PAV+, Resp Mech, Trending & LC. Unit is in good physical condition. Battery pack is working fine. The alarm display membrane is worn out.

Lot 592 Draeger 5400 Transport Incubator with Ferno Trolley, Air & O2 Hoses | SN ARNA

Serial No: ARNA-0012

Accessories: Ferno Trolley

Air & O2 hose

Condition Report: The incubator power up with the temperature control functional. The set temperature is maintained accordingly. Pneumatic system for O2 therapy and low suction were not tested. System working on mains, battery is not included. Ferno trolley is in good condition. System is in good overall condition.

Lot 593 Ethicon Mammutome Breast Biopsy System with Mammutome ST

Serial No: 2300118

Accessories: Mammutome ST Holster

Mammutome EX Biopsy Device (ultrasound)

Condition Report: System powers up with the LCD working. System requires all devices to be connected. Unable to do any further tests

Lot 594 Datex Ohmeda Aespire View Anaesthetic Machine with Monitor | DOM 2009 | SN

Manufactured: 2009

Accessories: Air, N2O & O2 hose

AGSS hose



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

Central Vacuum hose

Software: Version 6.20

Condition Report: The anaesthetic machine powers up and completes the initial system self test successfully. Pipeline pressure gauges, flowmeters and O2 flush are all functional. Anti-hypoxic feature and O2 basal flow are working fine. Aespire View Ventilator's LCD, rotary knob and key pad are working fine. All ventilation modes (VCV, PCV, SIMV/PSV, PSVPro, SIMV-PC & PCV-VG) are in good working order. The breathing assembly is in good condition. Flow sensors & O2 sensor are working fine. System was tested with medical gas and test lung, good ventilation cycles were delivered. The set tidal volume and breath frequency correspond to the resulting minute volume. The backup battery is working. The Aespire View is in good overall condition with minor scratches. One of the casters is not working properly. O2 & Air pin-index yoke is fitted at the back.

Lot 596 Olympus Stack Trolley, Sony Printer, Olympus CLV-260 Light Source, CV-260

Condition Report: The CV260 processor power up with indicators and function keys working LCD and keyboard are working fine. The processor was tested with a compatible videoscope and fair image was displayed on the 19" LCD monitor. The light source CLV 260 power with the xenon lamp and brightness control working fine. Sony printer UP-897MD, Panasonic DVD recorder LQ-MD800 power up ok but not fully tested. Sony digital capture unit failed to power up and the Sony colour video printer power up but with errors. Lot is in poor cosmetic condition, covered with dust and scuffs.

Lot 597 Haag Streit Bern 900 Slit Lamp with Eyepieces (Damaged)

Condition Report: The slit lamp power up with light beam intensity control and bulb working fine. Light beam patterns and filter are working fine. The optical system was tested with Test Rod, clear optic with good magnification and focus. The slit lamp mount roller is damaged on one side. Table height adjustment is working with few scratches on the base.

Lot 598 Toshiba Nemio 10 Ultrasound with Toshiba PVM-375AT Convex Array Transducer

Serial No: E2524090

Manufactured: 02 2002

Associated Equipment: Toshiba PVM-375AT Convex Array Transducer
Serial Number: TDA12Z2875
Manufactured: 12-2010

Mitsubishi P91 Printer

Sony SVO-9500MDP VCR

Condition Report: The ultrasound powers up and successfully completes its boot-up process. The CRT monitor, control panel function keys, trackball, TGC sliders and alphanumeric keyboard are in good working order. When tested using a phantom the probe produced



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

satisfactory image quality, resolution, grey-scale and accurate caliper measurement were obtained. The functional system modes are 2D and M-Mode. The probe has major crystals dropout on the edge and a small section towards the middle of the image. The transducer membrane is discoloured, the cable sheath has a small cut and the strain relief is loose on the side of the transducer connector. The printer and the foot switch are in good working order. The VHS recorder powers up and appears to be working but was not fully tested. The door release button of the printer is broken but it can be opened. One of the plastic cover of the transducer holder is missing. Overall the ultrasound is in satisfactory physical condition with minor signs of wear.

Lot 599 Sonosite iLook Portable Ultrasound with Transducer and Power Pack | DOM

Serial No: 038L5D

Manufactured: 05 2007

Accessories: Linear array probe

Software: ARN Ver: 10.80.100.025

Condition Report: Unit power up with mains with LCD and soft keys working fine. System function menus are accessible and working. Image displayed are of fair quality. Gain adjustments are working. The internal battery is low, was not working during test. System is in a fair overall condition.

Lot 705 Haag Streit 900 Slit Lamp | 4612

Serial No: 900.7.1.6347

Accessories: Eye piece set

Condition Report: The slit lamp power up with light beam intensity control and bulb working fine. Light beam patterns and filter are working fine. The optical system was tested with Test Rod, clear optic with very good magnification and focus. Table height adjustment is working with few scratches on the base.

Lot 709 Hill-Rom Airshields Resuscitaire | 11957

Serial No: L101757

Condition Report: The unit powers up and the displays, function keys are working fine. The overhead heater and light are working fine when tested in the pre-warm and manual modes. The baby mode was not fully tested, skin temperature probe is not included. The motorised height adjustment is in good working order. The pneumatics system include suction, Auxiliary flow & patient supply was not fully tested. Mechanical cot tilt and mattress are in good condition, with minor scratches on the handle. Lot in good and cosmetic physical condition. With O2 inlet only.

Lot 725 Olympus Colposcope | 04810

Serial No: 503020

Accessories: x2 Eye pieces GSWH20X-H-2/12,5

Condition Report: The light source power up with halogen lamp and brightness control in good working order. The fibre optic guide cable has been tapped at one end but the light transmission is good. The optical system is clear with good magnification and focus. One of the



UK Head Office: Gloucester Road Trading Estate Malmesbury Wiltshire SN16 9JT
+44 (0) 1666 822577 sales@hilditchgroup.com www.hilditchgroup.com

support arm lock handle is broken, otherwise the unit is in good physical and cosmetic condition.

Lot 762 Valley Lab Force 2 Diathermy with Foot Pedals

Serial No: F8K36258T

Manufactured: 08 1998

Associated Equipment: 1x Bipolar Footswitch
1x Monopolar Footswitches

Condition Report: Unit power up with the alphanumerical displays and function keys working fine. The display power up test completed. Patient return plate monitoring circuit is functional. The system was tested with a ESU analyser in both bipolar and monopolar functions and different modes. All the RF outputs were within the set limits. Unit includes a bipolar footswitch and a set of monopolar footswitches, which are both in good working condition. The unit is in fair physical condition with some minor scratches and marks, and slight uplifting of cover at the front.